

## High Power CO<sub>2</sub> Fractional Laser for Skin Resurfacing & Vaginal Rejuvenation

# FRA<sub>XIS</sub>

FRAXIS is a fractional CO<sub>2</sub> system that delivers precise and accurate energy to the targeted tissue with predictable ablation and less pain for optimal treatment result.

- **True 30W** High Power Tube
- **100 $\mu$ m** Ideal Beam Size
- **Multi shape** & Scanning Pattern
- **Less Pain** and Down Time

### INDICATIONS

#### Anti aging

- Tone & Texture
- Pore Reduction
- Wrinkle reduction
- Skin lifting & tightening
- Acne Scar
- Whitening
- Resurfacing

#### GynoLaser™

- Vaginal tightening
- Urinary Incontinence



# FRAxis

## High Power CO<sub>2</sub> Fractional Laser for Skin Resurfacing & Vaginal Rejuvenation

"Fraxis® CO<sub>2</sub> laser is a reference in terms of comfort and safety for treatment of various benign indications as well as profitability in aesthetic indications. The advantage of the Fraxis® CO<sub>2</sub> laser is its reliability of use, robustness, and power related to its RF emission tube."

- Dr. Michael Naouri, Dermatologist

"I have found the Fraxis very useful and essential for my clinic, and it has given a big boost to my practice. I would strongly recommend this device to all my colleagues."

- Dr. Shahid Javaid Akhtar, Dermatologist

Photos courtesy of Shahid Javaid Akhtar, M.D.



Before treatment



Six months after third Fraxis treatment

### SPECIFICATION

Laser Type	RF CO <sub>2</sub> Laser, All Metal Sealed Type
Laser Power	30W, CW/Pulse/Ultra mode: 0.1 ~ 30W
Laser Mode	10.6um(TEM00)
Pulse Duration	20 - 5,000us
Repetition	0.2 / 0.5 / 0.7 / 1s / 1.5s / 2sec / single
Overlap	1 - 10th
Distance	0.1 - 2.0mm
Treatment Area(mm)	1 x 1 - 20 x 20
Pixel Quantity	Up to 40, 401
Pixel Size	≥ 100 micron
Cooling	Air Cooling
Optical Guide	Articulated Arm
Aiming Beam	Semiconductor Diode laser/ Max 5mW/ 655nm

### OPTION

Surgical Handpiece 0.2 / 0.3 / 0.5 / 0.7 / 0.9 / 1.1 / 1.3mm



GYNOLASER

Non-surgical Laser treatment for vaginal tightening and urinary incontinence



▲ Disposable external slide

▲ 90° degree probe set



▲ Distance gauge



GYNOLASER



DISTRIBUIDOR OFICIAL

Avda. Dionisio Cifuentes, 34 · 33203 Gijón 985 30 12 34  
info@belium.es · www.beliummedical.es

**BELIUM**  
MEDICAL